

1       **Generation of anti-azoxystrobin monoclonal antibodies**  
2       **from regioisomeric haptens functionalized at selected**  
3       **sites and development of indirect competitive**  
4       **immunoassays**

5                               Javier Parra<sup>a</sup>, Josep V. Mercader<sup>a</sup>, Consuelo Agulló<sup>b</sup>,

6                               Antonio Abad-Somovilla,<sup>b,\*</sup> and Antonio Abad-Fuentes<sup>a,\*</sup>

7  
8       <sup>a</sup> *Department of Biotechnology, IATA–CSIC, Agustí Escardino 7, 46980 Paterna, València, Spain*

9       <sup>b</sup> *Department of Organic Chemistry, Universitat de València, Doctor Moliner 50, 46100 Burjassot,*  
10       *València, Spain*

11  
12       \* Corresponding authors: Dr. A. Abad-Fuentes, Tel. +34-963900022; Fax +34-963636301; e-mail  
13       [aabad@iata.csic.es](mailto:aabad@iata.csic.es) and Dr. A. Abad-Somovilla, Tel. +34-963544509; Fax +34-963544328; e-mail  
14       [antonio.abad@uv.es](mailto:antonio.abad@uv.es)

15

16 **Abstract**

17 Azoxystrobin is a modern strobilurin fungicide used around the world to combat prime diseases  
18 affecting highly valuable crops. Accordingly, residues of this chemical are frequently found in food,  
19 even though mostly under maximum tolerated levels. We herein describe the development of an indirect  
20 competitive immunoassay for the determination of azoxystrobin residues. A panel of monoclonal  
21 antibodies displaying subnanomolar affinity to azoxystrobin was generated using, as immunizing  
22 haptens in mice, four functionalized derivatives carrying the same spacer arm located at different  
23 rationally-chosen positions. This collection of antibodies was thoroughly characterized with  
24 homologous and heterologous antigens, and the immunoassay consisting of monoclonal antibody  
25 AZo6#49 and the coating conjugate OVA–AZb6, which displayed an  $IC_{50}$  value of  $0.102 \mu\text{g L}^{-1}$  and a  
26 LOD of  $0.017 \mu\text{g L}^{-1}$ , was eventually optimized. The response to different pH and ionic strength  
27 conditions of the specific assay was studied using a biparametric approach. In addition, the influence of  
28 Tween 20 and organic solvents over the assay parameters was also evaluated. After optimization, the  
29 developed immunochemical assay was applied to the analysis of azoxystrobin in spiked juices of  
30 relevant fruits and vegetables, showing excellent recoveries between 2 and  $500 \mu\text{g L}^{-1}$ .

31

32

33 **Keywords**

34 competitive ELISA; site heterology; Landsteiner's principle; fungicide; strobilurins; food safety

35

36

37 **Abbreviations**

38 AZ, azoxystrobin; BSA, bovine serum albumin; CR, cross-reactivity; DMEM, Dulbecco's modified  
39 Eagle's medium; ELISA, enzyme-linked immunosorbent assay; FBS, fetal bovine serum; HAT,  
40 hypoxanthine–aminopterin–thymidine; HFCS, hybridoma fusion and cloning supplement; HT,  
41 hypoxanthine–thymidine; i-cELISA; indirect competitive ELISA; LOD, limit of detection; LOQ, limit of  
42 quantification; mAb, monoclonal antibody; MRL, maximum residue limit; OVA, ovalbumin; PBS,  
43 phosphate-buffered saline; PBST, phosphate-buffered saline with Tween 20; RAM–HRP, polyclonal  
44 rabbit anti-mouse immunoglobulin peroxidase conjugate.

45

## 46 1. Introduction

47 Strobilurins comprise a family of powerful fungicides. Their discovery was inspired by the  
48 identification of a group of active natural products displaying a potent antifungal activity in the  
49 mushrooms *Strobilurus tenacellus* and *Oudemansiella mucida* [1,2]. These agrochemicals exert their  
50 fungicide action through binding to the Q<sub>o</sub> site of cytochrome b in the mitochondrial electron transport  
51 chain [3]. Although strobilurins show much lower toxicity to humans and non-target organisms than the  
52 traditional anti-mycotic substances, some of them are toxic to aquatic organisms [4]. Since their  
53 commercialization in the late 1990s, the strobilurin family has been steadily growing. Currently, 17  
54 active compounds have been developed, although just 9 of them have received an ISO name and are  
55 therefore of global commercial relevance [5,6]. One of the two first strobilurins to be patented was  
56 azoxystrobin (AZ, Fig. 1), which was soon after approved worldwide for disease control in most cereals,  
57 fruits, and vegetables, being currently registered for use in nearly 120 crops in around 100 countries. At  
58 present, this active principle is being massively sold in a variety of formulations under different trade  
59 names. As a matter of fact, nearly 4000 tones of AZ was used worldwide in 2009, with global annual  
60 sales over \$1 billion, which makes it the world's leading proprietary fungicide [7].

61 Nowadays, the monitoring of pesticide residues is compulsory for governments and private  
62 corporations involved in food processing. The European maximum residue limits (MRLs) for AZ in  
63 most foodstuffs range between 0.05 and 5 mg kg<sup>-1</sup> [8]. According to the European Pesticide Monitoring  
64 Programs [9], AZ was among the most frequently found pesticides – around 5% of the food samples  
65 being analyzed contained residues at or below the MRL. Residues of this fungicide were mainly found  
66 in grapes, strawberries, peaches, tomatoes, leeks, head cabbages, peppers, peas, and lettuces. In  
67 addition, AZ residues have been found in locally purchased mangoes at levels between 12.7 and 55.8  
68 µg kg<sup>-1</sup> (16 samples, 100%) [10], and in commercial white and red wines at levels between 0.2 and 2.9  
69 µg L<sup>-1</sup> (11 samples, 82%) [11]. The analysis of AZ in foods is usually carried out by gas  
70 chromatography coupled to mass spectrometry detection (GC–MS) [12–14]. Other approaches have also

71 been published using fluorescence detection [15], high-performance liquid chromatography with diode  
72 array detection [16], desorption electrospray ionization mass spectrometry [17], ultrahigh-performance  
73 liquid chromatography tandem mass spectrometry [18], direct analysis in real time coupled to time-of-  
74 flight mass spectrometry [19], or liquid chromatography with electrospray tandem mass spectrometry  
75 [20].

76 For many applications, antibody-based bioanalytical techniques represent an attractive alternative to  
77 classical chromatographic methods. Immunoassays are generally considered as easy-to-use, affordable,  
78 rapid, sensitive, specific, and environmentally friendly analytical approaches, especially indicated for  
79 laboratories involved in environmental or food quality and safety programs [21–23]. Furthermore, due  
80 to their high versatility, immunochemical methods can easily be adapted to different particular  
81 analytical needs. However, immunoassay development demands raising suitable antibodies against the  
82 target analyte and evaluation of the produced immunoreagents into appropriate assay formats. The  
83 successful generation of specific and sensitive antibodies against small organic molecules, such as AZ,  
84 is greatly dependent upon the proper design of the immunizing and assay haptens, which must be  
85 coupled to a carrier protein. Following Landsteiner’s seminal works [24], the idea that any modification  
86 of the target molecule should be introduced at a distal site from unique determinant moieties has been  
87 accepted as a general rule to prepare immunogens. However, this so-called “distal position” is  
88 sometimes uncertain.

89 In 2006, Furzer et al. [25] described the production of polyclonal antibodies against AZ using the  
90 acidic form of the molecule for direct conjugation. However, monoclonal antibodies (mAbs) are more  
91 attractive immunoreagents to industry because they constitute a defined reagent and a constant supply  
92 can be guaranteed. Up to now, we have reported the synthesis of haptens, the production of mAbs, and  
93 the development of immunoassays for different strobilurin pesticides, *i.e.* kresoxim-methyl,  
94 trifloxystrobin, pyraclostrobin, and picoxystrobin [26–29], and also a comprehensive strategy for the  
95 synthesis of four derivatives of AZ functionalized at different positions with the same spacer arm has  
96 recently been reported by our group [30]. For the present study, mice were immunized with

97 bioconjugates of those regioisomeric synthetic haptens, and a collection of mAbs against AZ was  
98 generated. These mAbs were exhaustively characterized and an enzyme-linked immunosorbent assay  
99 (ELISA) was developed and optimized in the conjugate-coated indirect competitive format (i-cELISA).  
100 Finally, evaluation of the selected immunoassay was undertaken by determining AZ recoveries in  
101 fortified juice samples of relevant commodities.

## 102 **2. Materials and methods**

### 103 *2.1. Chemicals and instrumentation*

104 Analytical-grade AZ (methyl (E)-2-{2-[6-(2-cyanophenoxy)pyrimidin-4-yloxy]phenyl}-3-  
105 methoxyacrylate) (CAS number 131860-33-8, MW 403.4 g mol<sup>-1</sup>) was kindly provided by Syngenta  
106 AG (Basel, Switzerland). Sepharose HiTrap Protein G HP columns used for antibody purification were  
107 purchased from General Electric Healthcare (Uppsala, Sweden). Polyclonal rabbit anti-mouse  
108 immunoglobulin peroxidase conjugate (RAM-HRP) was from Dako (Glostrup, Denmark). The  
109 immunoglobulin isotype was determined using the Mouse MonoAb-ID kit from Invitrogen (Carlsbad,  
110 CA, USA). Freund's adjuvants, thimerosal (sodium ethylmercurithiosalicylate), ovalbumin (OVA), and  
111 *o*-phenylenediamine were purchased from Sigma-Aldrich (Madrid, Spain). Bovine serum albumin  
112 (BSA) and Hybridoma Fusion and Cloning Supplement (HFCS) were from Roche Applied Science  
113 (Mannheim, Germany). P3-X63-Ag8.653 mouse plasmacytoma cell line was from the European  
114 Collection of Cell Cultures (Wiltshire, UK). Gentamicine, hipoxanthine-thimidine (HT) and  
115 hipoxanthine-aminopterin-thimidine (HAT) solutions were obtained from Gibco BRL (Paisley,  
116 Scotland). Cell culture media (high-glucose Dulbecco's Modified Eagle's Medium, DMEM),  
117 polyethylene glycol Hybri-Max (PEG 1500), fetal bovine serum (FBS), 200 mM alanyl-glutamine  
118 solution, red blood cell lysing buffer Hybri-Max, and MEM non-essential amino acid solution were  
119 from Sigma-Aldrich (Madrid, Spain). Culture plasticware and Costar flat-bottom high-binding  
120 polystyrene ELISA plates were from Corning (Corning, NY, USA). ELISA absorbances were read with

121 a PowerWave HT microplate reader and microplates were washed with an ELx405 microplate washer,  
122 both from BioTek Instruments (Winooski, VT, USA).

123 The synthesis of the four functionalized derivatives of AZ used in this work (haptens AZa6, AZb6,  
124 AZc6, and AZo6, Fig. 1), as well as the preparation of the immunizing and coating conjugates, has been  
125 previously described [30].

## 126 2.2. Monoclonal antibody production

127 Animal manipulation was performed in compliance with the Spanish laws and guidelines (RD  
128 1201/2005 and law 32/2007) and according to the European Directive 2003/65/EC concerning the  
129 protection of animals used for experimental and other scientific purposes. Four groups of four  
130 BALB/cByJ female mice (8–10 weeks old) were immunized by intraperitoneal injections with BSA–  
131 AZa6, BSA–AZb6, BSA–AZc6, or BSA–AZo6 conjugates. Doses consisted of 100 µg of protein  
132 conjugate in Freund's adjuvant (complete for the first injection and incomplete for the second and the  
133 third ones). The antiserum from each mouse was obtained by submandibular bleeding 9–10 days after  
134 the third injection. Sera were diluted 1/5 with PBS (10 mM sodium phosphate buffer, 140 mM NaCl,  
135 pH 7.4) containing 0.01% (w/v) thimerosal and stored at 4 °C in amber glass vials. The antiserum titer  
136 was calculated as the required dilution to obtain an absorbance of 1.0 under standard ELISA conditions  
137 using the homologous conjugate (a conjugate with the same hapten that was used for immunization) at  
138 1.0 µg mL<sup>-1</sup>. After a resting period of at least 3 weeks from the last injection with adjuvant and four  
139 days before cell fusion, mice received a booster intraperitoneal injection of 100 µg of protein conjugate  
140 in 200 µL of PBS.

### 141 2.2.1. Cell fusion and culture

142 P3-X63/Ag 8.653 murine myeloma cells were cultured in high-glucose DMEM supplemented with  
143 2 mM alanyl-glutamine, 1 mM MEM non-essential amino acids, and 25 µg mL<sup>-1</sup> gentamicin (referred  
144 to as s-DMEM) and containing 10% (v/v) FBS. Mouse spleen lymphocytes were fused with myeloma

145 cells at a 4:1 ratio using 1 mL of PEG 1500 as the fusing agent. The fused cells were distributed in 96-  
146 well culture plates at a density of  $1.5 \times 10^5$  lymphocytes per well in 100  $\mu$ L s-DMEM with 15% FBS.  
147 Twenty-four hours after plating, 100  $\mu$ L of selection medium (s-DMEM supplemented with HAT) with  
148 20% FBS and 1% (v/v) HFCS was added to each well.

### 149 *2.2.2. Hybridoma selection and cloning*

150 Twelve days after fusion, hybridoma culture supernatants were assayed following a sequential double-  
151 screening process. First, each culture supernatant was analyzed in parallel with and without competitor  
152 by i-cELISA as described by Abad et al. [31]. The concentration of AZ in solution was 500 nM in the  
153 screening of the first cell fusion of each set of mice, and 100 nM in successive fusion experiments. The  
154 signal in non-competitive conditions (absence of analyte) was compared with the competitive one when  
155 AZ was used as competitor, and the ratio of both absorbances was used as the criterion to identify good  
156 antibody-secreting clones. In addition, competitive ELISAs were performed using serial dilutions of the  
157 culture supernatant from those wells that afforded saturated signals in the first screening experiment.  
158 The selective pressure was increased in the second assay using lower AZ concentrations, and ELISA  
159 plates were coated with antigens at  $0.1 \mu\text{g mL}^{-1}$  to favor competition. The selected hybridomas were  
160 cloned by limiting dilution in cloning medium (s-DMEM containing 20% FBS and supplemented with  
161 HT and 1% HFCS). Stable antibody-producing clones were expanded and cryopreserved in liquid  
162 nitrogen.

### 163 *2.2.3. Purification of monoclonal antibodies*

164 Immunoglobulins were purified from late stationary-phase culture supernatants by ammonium sulfate  
165 precipitation and protein G affinity chromatography following column manufacturer's instructions. A  
166 fraction of the purified mAb was diluted 1:1 with PBS containing 1% (w/v) BSA and 0.01% (w/v)  
167 thimerosal and stored at 4 °C in amber glass vials for daily usage. The remaining mAb solution was  
168 stored at 4 °C as ammonium sulfate precipitate.



## 170 2.3. Enzyme-linked immunosorbent assays

171 Ninety-six-well polystyrene ELISA plates were coated with 100  $\mu\text{L}$  per well of OVA conjugate  
172 solution in 50 mM carbonate–bicarbonate buffer, pH 9.6, by overnight incubation at room temperature.  
173 Next day, coated plates were washed four times with 0.15 M NaCl containing 0.05% (v/v) Tween 20,  
174 and received 50  $\mu\text{L}$  per well of AZ standard solutions plus 50  $\mu\text{L}$  per well of antibody diluted in PBST  
175 (PBS containing 0.05% (v/v) Tween 20). AZ standard curves were prepared by serial dilution in PBS  
176 using borosilicate glass tubes from a concentrated AZ stock solution in anhydrous *N,N*-  
177 dimethylformamide. A blank was included in each curve. The immunological reaction took place during  
178 1 h at room temperature. After washing as above, plates received 100  $\mu\text{L}$  per well of a 1/2000 dilution  
179 of RAM–HRP conjugate in PBST, and they were incubated 1 h at room temperature. Plates were  
180 washed again and the signal was generated by adding 100  $\mu\text{L}$  per well of freshly prepared enzyme  
181 substrate solution (2 mg  $\text{mL}^{-1}$  of *o*-phenylenediamine and 0.012% (v/v)  $\text{H}_2\text{O}_2$  in 25 mM sodium citrate  
182 and 62 mM sodium phosphate buffer, pH 5.4). The enzymatic reaction was stopped after 10 min at room  
183 temperature by addition of 100  $\mu\text{L}$  per well of 2.5 M sulfuric acid. The absorbance was immediately  
184 read at 492 nm using 650 nm as reference wavelength.

185 Competitive curves were obtained by plotting mean absorbance values versus the logarithm of analyte  
186 concentration. The resulting experimental values were fitted to a four-parameter logistic equation using  
187 the SigmaPlot software package from SPSS Inc. (Chicago, IL, USA). From those curves, assay  
188 sensitivity was estimated as the concentration of analyte at the inflection point of the fitted curve,  
189 typically corresponding to a 50% inhibition ( $\text{IC}_{50}$ ) of the maximum absorbance ( $A_{\text{max}}$ ) if the background  
190 signal approaches to zero. The limit of detection (LOD) was estimated as the concentration of AZ that  
191 provided a 10% inhibition of the  $A_{\text{max}}$  ( $\text{IC}_{10}$ ), whereas the limit of quantification (LOQ) was established  
192 at the  $\text{IC}_{20}$  and the working range between the  $\text{IC}_{20}$  and  $\text{IC}_{80}$  values. Cross-reactivity (CR) was  
193 calculated according to the formula:

194 
$$CR = [IC_{50} (AZ) / IC_{50} (\text{compound})] \times 100$$

195 *2.4. Buffer ionic strength and pH studies*

196 Biparametric analysis of assay response to diverse pH and ionic strength conditions was performed  
197 using statistical software from Minitab Inc (State College, PA, USA). Axial points were fixed at 5.5 and  
198 9.5 for pH, and at 50 and 300 mM for the ionic strength. The number and composition of buffers  
199 indicated by the software were prepared as follows. A 40 mM trisodium citrate, 40 mM disodium  
200 hydrogen phosphate, 40 mM Tris solution was made (pH 9.9,  $I = 360$  mM), and different volumes of  
201 5 M HCl and 2 M NaCl were added to aliquots of that solution in order to achieve buffers with different  
202 pH and ionic strength values. The ionic strength of each solution was determined according to the  
203 formula

204 
$$I = \frac{1}{2} \sum c_i z_i^2$$

205 where  $I$  is the ionic strength,  $c$  is the concentration of each ion at equilibrium, and  $z$  is its charge. After  
206 pH and  $I$  adjustment at each desired value, Tween 20 and deionized water were added, so the final  
207 composition of every buffer was 20 mM citrate, 20 mM phosphate, and 20 mM Tris containing 0.05%  
208 Tween 20 and different amounts of NaCl. Competitive assays were performed with AZ standards in  
209 deionized water.

210 *2.5. Sample analysis*

211 Commercial juices from tomatoes, grapes, peaches, and bananas were acquired from a local  
212 supermarket. Analysis of those samples by GC-MS (see the [Supplementary Data File](#) for  
213 chromatographic conditions) proved the absence of AZ residues at levels higher than  $20 \mu\text{g kg}^{-1}$  (LOD  
214 of the employed experimental procedure). Competitive assays were performed in plates coated with a  
215 solution of OVA-AZb6 at  $1 \mu\text{g mL}^{-1}$ . Each coated well received 50  $\mu\text{L}$  of AZ standard solutions in PBS  
216 or 50  $\mu\text{L}$  of sample fortified with AZ and diluted in PBS (1/25, 1/100, and 1/500), plus 50  $\mu\text{L}$  of a 60 ng  
217  $\text{mL}^{-1}$  solution of mAb AZo6#49 prepared in PBST. Each sample was analyzed in triplicate wells, and

218 the mean absorbance values were interpolated into the standard curve run in the same plate, also in  
219 triplicate wells. Unfortified samples were also diluted in PBS as mentioned above and included as  
220 controls.

### 221 **3. Results and discussion**

#### 222 *3.1. Immunological response*

223 In order to elicit antibodies suitable for the development of sensitive immunochemical assays for  
224 small organic molecules, the immunizing haptens must mimic as much as possible the structure and  
225 electronic distribution of the parent analyte. AZ is a very flexible compound, so it is rather uncertain to  
226 predict the linker position that would limit to a lesser extent the conformational freedom of the  
227 molecular framework, and simultaneously afford a maximum exposure of the target to the immune  
228 system. In this study, four rationally-designed haptens, namely AZa6, AZb6, AZc6, and AZo6 (Fig. 1),  
229 which incorporated the same saturated hydrocarbon spacer arm at different positions of the skeleton  
230 through C–C single bonds, were evaluated for the production of mAbs. The syntheses of those  
231 regioisomeric functionalized haptens of AZ and their conjugation to carrier proteins have been  
232 described elsewhere [30].

233 The immunization process was confirmed and the performance of the immunogens was assessed  
234 using antisera collected after the third boost. At this point, all of the animals showed an outstanding  
235 response, with antisera showing high titers ( $\sim 10^5$ ) when they were evaluated by ELISA using the  
236 homologous coating conjugate. In addition, all antisera bound all heterologous coating antigens  
237 (conjugates bearing a hapten different to that of the immunogen) with good titers despite the diversity in  
238 spacer arm attachment sites, even at low conjugate concentrations ( $0.1 \mu\text{g mL}^{-1}$ ). Concerning affinity,  
239 the obtained  $\text{IC}_{50}$  values for AZ were also remarkably low for mouse antisera, particularly those of mice  
240 immunized with haptens AZa6 and AZb6. Moreover, some heterologous combinations resulted in up to  
241 6-fold reduction of the  $\text{IC}_{50}$  values. As listed in Table S1 (Supplementary Data File), the lowest  $\text{IC}_{50}$

242 values were achieved with AZa6-type antisera using OVA–AZb6 for coating. In general, OVA–AZb6  
243 conjugate was the best heterologous coating antigen for all antisera, probably because of the central  
244 position of the spacer arm. Anyhow, the four evaluated haptens gave rise to an excellent response in  
245 mice, so all of them were deemed appropriate candidates to address hybridoma technology.

### 246 *3.2. Generation of monoclonal antibodies*

247 At least three cell fusions were performed with splenocytes from mice immunized with each AZ  
248 derivative. A sequential double-screening procedure was carried out after the cell fusion for hybridoma  
249 selection, as described in the section 2.2.2. In the first screening, two simultaneous assays – one with  
250 and one without free analyte – were performed, which helped to identify those antibody-producing  
251 clones that bound not only the coating conjugate but also the free analyte. The observed high number of  
252 positive clones (wells containing antibodies that recognized the coating conjugate) that were usually  
253 found after each cell fusion could already be envisaged from the high titers exhibited by mouse antisera.  
254 In most cases, also a high number of competitive clones (wells containing antibodies that recognized the  
255 free analyte) were retrieved in this first screening assay (Table S2). With the aim of efficiently ranking  
256 hybridomas for further cloning and expansion according to the affinity of the secreted antibodies, a  
257 second screening was routinely carried out following a competitive checkerboard assay. Remarkably,  
258 high-affinity antibody-producing hybridomas were rather homogeneously derived from all of the  
259 immunizing haptens, in keeping with the results observed in competitive experiments with the mouse  
260 antisera. In total, 37 hybridoma cell lines were cloned and stabilized; they all produced antibodies of the  
261 IgG<sub>1</sub> isotype with  $\kappa$  light chains, with the only exception of mAb AZb6#43 (IgG<sub>2a</sub>,  $\kappa$ ).

### 262 *3.3. Evaluation of antibody affinity*

263 A straightforward strategy for immunoreagent assessment and affinity determination was followed.  
264 Each mAb, at six different concentrations (30–1000 ng mL<sup>-1</sup>), was tested against the four available  
265 coating conjugates at two concentrations (0.1 and 1.0  $\mu$ g mL<sup>-1</sup>) under competitive conditions using

266 8-point AZ standard curves ( $10^{-4}$ – $10^2$  nM), including a blank. All combinations were tested at least  
267 three times. This study provided, for each pair of immunoreagents, a set of 12 inhibition curves  
268 displaying different  $A_{\max}$  values, which allowed a well-founded selection of the best combinations of  
269 mAb and assay conjugate for further ELISA development. Concerning antigen recognition, all of the  
270 mAbs were able to bind the homologous conjugate with high avidity. In fact, coating antigens had to be  
271 used at or below  $0.1 \mu\text{g mL}^{-1}$  in homologous combinations to afford optimum competitive curves. With  
272 regard to heterologous coating conjugates, unexpected results were found in some cases. Thus, most  
273 AZa6-type mAbs were unable to bind OVA–AZc6, whereas the opposite was not true, *i.e.* AZc6-  
274 derived mAbs recognized pretty well OVA–AZa6 as coating antigen. A similar finding was observed  
275 with haptens AZb6 and AZa6; whereas most AZb6-type mAbs did not recognize OVA–AZa6, 7 out of  
276 11 AZa6-type mAbs did bind OVA–AZb6. Finally, we found that conjugate OVA–AZo6 was strongly  
277 recognized by all of the antibodies, even by those generated from immunizing haptens with the linker at  
278 a very different site. Most likely, the many degrees of freedom of the toxophore moiety in AZo6 could  
279 probably explain its universal character as coating antigen.

280 Table 1 summarizes the antibody and coating conjugate concentrations, as well as the curve  
281 parameters, for a short selection of mAbs at those particular immunoreagent concentrations that  
282 afforded the lowest  $\text{IC}_{50}$  with an  $A_{\max}$  value between 0.8 and 1.5 (see Table S3 for the rest of mAbs).  
283 Overall, antibodies displaying outstanding affinity to AZ ( $\text{IC}_{50}$  values in the subnanomolar range) were  
284 obtained from any of the four immunizing conjugates. The antigen OVA–AZb6, with the spacer arm  
285 located at the central ring of the AZ skeleton, emerged as a very convenient heterologous coating  
286 conjugate for the development of sensitive assays, as it was already observed with mouse antisera.  
287 Slopes were an issue of concern from this study because very steep inhibition curves were observed  
288 with many immunoreagent combinations, which would result in assays with rather narrow working  
289 ranges. Accordingly, this curve parameter received special attention in the further selection of  
290 immunoreagent combinations.

### 291 3.4. Assessment of antibody specificity

292 The binding properties of each mAb towards relevant members of the strobilurin family (kresoxim-  
293 methyl, trifloxystrobin, picoxystrobin, pyraclostrobin, dimoxystrobin, oryastrobin, and  
294 metominostrobin) were studied using the homologous OVA conjugate as coating antigen. No inhibition  
295 was observed up to 1  $\mu\text{M}$  in any case, so all mAbs, irrespective of their origin, were highly specific to  
296 AZ. This finding was not actually unexpected because strobilurin fungicides, although belonging to the  
297 same agrochemical group, barely share a minimum part of the whole structure (Fig. S1), and it is well-  
298 known that mAbs are biomolecules most often displaying exquisite specificity to their target, being able  
299 to discriminate between chiral molecules and geometric isomers [32,33].

300 The fungicide activity of AZ strongly relies on the *E*-isomer, which in fact makes up for more than  
301 98% of the technical product, whereas the *Z*-isomer, which is mainly generated by photochemical  
302 transformation, accounts for up to 25% of the observed degradation products [34,35]. During hapten  
303 synthesis, especial precautions were taken to preserve the *E* conformation of the toxophore in all  
304 derivatives. Accordingly, we thought it was worthwhile to challenge our large collection of mAbs with  
305 the *Z*-isomer in order to further explore the fine specificity of their binding sites and to study the  
306 relationship between their stereoselectivity and the derivatization position in the parental hapten. The  
307 *Z*-isomer of AZ was prepared as described by Clough et al. [36], and competitive curves with both  
308 isomers were simultaneously run in the homologous format with the whole collection of mAbs. The  
309 most stereoselective mAbs were those coming from immunizing haptens AZo6 and AZb6, showing  
310 mean CR values for each group of mAbs of 3.8% ( $n = 12$ ) and 3.5% ( $n = 8$ ), respectively (Fig. 2). On  
311 the contrary, mAbs derived from the immunizing hapten AZa6 were less able to discriminate between  
312 both isomers, with a mean CR value for the *Z*-isomer of 40.3% ( $n = 11$ ). Interestingly, this group  
313 included two of the antibodies that displayed extreme and opposite selectivity patterns, *i.e.* mAb  
314 AZa6#21, whose CR value for the *Z*-isomer was 1.0%, and mAb AZa6#210, with a CR value of 96.7%  
315 (Fig. S2). Regarding antibodies that came from the immunizing hapten AZc6, a mean CR value for the

316 Z-isomer of 19.0% was found. In brief, stereoselective antibodies (CR < 2%) were generated from all  
317 immunizing haptens, but thanks to the large set of available mAbs, we could observe that haptens AZo6  
318 and AZb6 were more adequate for this purpose. Otherwise, hapten AZa6 emerged as a very well suited  
319 derivative for the production of mAbs displaying equivalent affinities for both geometric isomers.

### 320 3.5. ELISA development and characterization

321 Considering that all of the produced antibodies displayed high affinity and specificity to AZ, some  
322 immunoreagents were reevaluated in order to ensure the selection of the best assays. Thus, those  
323 combinations that provided IC<sub>50</sub> values below 0.5 nM and slopes higher than -1.4 in the initial  
324 evaluation were tested again using adjusted mAb and antigen concentrations to afford A<sub>max</sub> values close  
325 to 1.0. Table 2 lists the main curve parameters that were reached with the best antibody that was derived  
326 from each type of immunogen. The most sensitive immunoassays were achieved with mAbs AZa6#11  
327 and AZo6#49, using in both cases OVA–AZb6 as coating conjugate. Taking into account IC<sub>50</sub> values,  
328 slopes, robustness, and cell growth and antibody-producing capacity in culture of the parental  
329 hybridoma cell line, mAb AZo6#49 together with conjugate OVA–AZb6 were selected for further assay  
330 characterization.

331 The influence of buffer pH and ionic strength was investigated following a biparametric approach. AZ  
332 concentrations from 100 nM to 1 pM were prepared in deionized water and the mAb was diluted in 20  
333 mM buffers of varying pH and *I* values. Fig. 3 shows the variation of the A<sub>max</sub> and IC<sub>50</sub> values of the  
334 resulting inhibition curves with regard to these two physicochemical parameters. High pH and *I* values  
335 decreased the A<sub>max</sub>, but they did not significantly modify the IC<sub>50</sub> inside the assayed range, even though  
336 slightly better IC<sub>50</sub> values were obtained in buffers with extreme pH and ionic strength. Therefore, the  
337 best compromise was to keep the ELISA conditions at pH and *I* values equivalent to those of the PBS  
338 (pH 7.4 and *I* = 162 mM at 25 °C). The moderate dependence of the assay IC<sub>50</sub> over the pH and salt  
339 concentration, together with the particularly high sensitivity of the described ELISA, should permit the

340 neutralization of the acidic character of most fruits and vegetables by simple dilution with a regular  
341 assay buffer.

342 Tween 20 is a common additive used in ELISA to avoid unspecific interactions between proteins and  
343 with the solid support. However, it is known that this detergent may exert a negative influence over the  
344 signal and sensitivity of the immunoassay [26]. As a further step in the optimization of the selected i-  
345 cELISA for AZ, competitive assays were run in PBS containing different amounts of Tween 20. From  
346 this study, it was observed that the usual detergent concentration (0.05%) was optimum for this assay  
347 (Fig. S3). Finally, the effect of some commonly employed organic solvents in pesticide residue  
348 extraction over the parameters of the inhibition curve was evaluated (Fig. S4). Solvent concentrations at  
349 or above 10% severely compromised the performance of the assay, in particular with acetonitrile and  
350 acetone. The best tolerated solvent was methanol; concentrations below 5% did not significantly change  
351 the curve parameters, and even acceptable  $A_{\max}$  and  $IC_{50}$  values were obtained with 10% methanol.  
352 Table 3 shows the main analytical features of the optimized competitive ELISA.

### 353 3.6. Food analysis

354 AZ is extensively used to fight botrytis in high-valuable crops, including tomatoes, peaches, grapes,  
355 and bananas, so the performance of the developed i-cELISA was evaluated in juices from those food  
356 commodities. Firstly, a study of the matrix effects over the developed immunoassay was undertaken by  
357 running AZ standard curves prepared in juice diluted in PBS. It was observed that a direct 1/50 dilution  
358 of the juice in assay buffer practically reduced the interferences down to negligible levels (Fig. 4). Next,  
359 the four juices were fortified with varying amounts of AZ and measured with the optimized i-cELISA  
360 after a simple dilution in buffer; three dilution factors were tested to better evaluate assay performance.  
361 As shown in Table 4, good to excellent recoveries over more than two orders of magnitude of analyte  
362 concentration could be obtained even from just a single dilution (factor 1/100). In the four analyzed  
363 commodities,  $2 \mu\text{g L}^{-1}$  of AZ could be adequately measured by direct dilution in PBS, which is a  
364 concentration 10 times below the lowest MRL that has been laid down in the international legislation.



365 Actually, the newly developed ELISA for AZ is currently being applied as a screening procedure in real  
366 samples within the framework of an ongoing project aimed at analyzing the presence of fungicide  
367 residues in market commodities by immunochemical methods. So far, two strawberry samples out of 16  
368 bought at local groceries have been found positive by the developed immunoassay. One of those  
369 samples contained AZ residues at a concentration of 620  $\mu\text{g kg}^{-1}$ . This sample was further analyzed by  
370 GC–MS, which provided a value of 630  $\mu\text{g kg}^{-1}$  (see Fig. S5 in the [Supplementary Data File](#) for the  
371 GC–MS chromatogram of that sample). The other positive sample contained traces of AZ according to  
372 both ELISA and the reference method (15–20  $\mu\text{g kg}^{-1}$ ).

#### 373 **4. Conclusions**

374 A collection of mAbs to AZ was generated from four haptens functionalized with the same linker  
375 length at different sites. Interestingly, antibodies showing high affinity to the target analyte were  
376 commonly found independently of the derivatization site in the immunizing hapten. However, a  
377 dependence of the antibody stereoselectivity upon the tethering site of the immunogen was envisaged. A  
378 detailed study of the antibody capacity to bind site-heterologous haptens revealed that OVA–AZo6 was  
379 recognized by all of the antibodies. Also, AZo6-derived immunoglobulins were the group of antibodies  
380 that better bound heterologous conjugates. In addition, the OVA–AZb6 conjugate emerged as the most  
381 appropriate antigen for the development of heterologous assays. Little influence of pH, Tween 20, and  
382 salt concentration was noticed over the optimized assay using mAb AZo6#49 and OVA–AZb6. Finally,  
383 the developed immunoassay showed an excellent performance in recovery studies with spiked juices,  
384 allowing the simple and rapid analysis of numerous samples in just an 8-hour working day. Particularly,  
385 this AZ ELISA could determine the target fungicide in fruit juices at levels as low as 2  $\mu\text{g L}^{-1}$ , which  
386 satisfactorily compares with instrumental chromatographic techniques.

387

388

389 **Acknowledgments**

390 We thank Laura López Sánchez and Ana Izquierdo Gil for excellent technical assistance.

391 This work was supported by *Ministerio de Educación y Ciencia* (AGL2006-12750-C02-01/02/ALI)  
392 and cofinanced by FEDER funds. J.P. and J.V.M. were hired by the CSIC, the former under a  
393 predoctoral *I3P* contract and the latter under a *Ramón y Cajal* postdoctoral contract, both of them  
394 financed by the Spanish *Ministerio de Ciencia e Innovación* and the European Social Fund.

395 Limited amounts of the immunoreagents described in this paper are available upon request for  
396 evaluation.

397 **Appendix A. Supplementary data**

398 Supplementary material associated with this article can be found, in the online version, at doi:

399

400 **References**

- 401 [1] T. Anke, F. Oberwinkler, W. Steglich, G. Schramm, *J. Antibiot.* 30 (1977) 806–810.
- 402 [2] H. Sauter, W. Steglich, T. Anke, *Angew. Chem. Int. Ed.* 38 (1999) 1328–1349.
- 403 [3] H. Balba, *J. Environ. Sci. Health B* 42 (2007) 441–451.
- 404 [4] D.M. Bartlett, J.M. Clough, J.R. Godwin, A.A. Hall, M. Hamer, B. Parr-Dobrzanski, *Pest. Manag. Sci.* 58  
405 (2002) 649–662.
- 406 [5] W. Huang, P.L. Zhao, C.L. Liu, Q. Chen, Z.M. Liu, G.F. Yang, *J. Agric. Food Chem.* 55 (2007) 3004–  
407 3010.
- 408 [6] Compendium of pesticide common names at <http://www.alanwood.net/pesticides/index.html>. Accessed  
409 April 14th, 2011.
- 410 [7] J. Atkin, UBS Basic Materials Conference (2010) at [http://www2.syngenta.com/en/investor\\_relations/  
411 event\\_calendar.html](http://www2.syngenta.com/en/investor_relations/event_calendar.html). Accessed April 14th, 2011.
- 412 [8] European Union Pesticide Residues Database at [http://ec.europa.eu/sanco\\_pesticides/public/index.cfm](http://ec.europa.eu/sanco_pesticides/public/index.cfm).  
413 Accessed April 14th, 2011.
- 414 [9] European Food Safety Authority at <http://www.efsa.europa.eu/en/scdocs/scdoc/305r.htm>. Accessed April  
415 14th, 2011.
- 416 [10] A. Menezes, F.N. dos Santos, P.A.D. Pereira, *Talanta* 81 (2010) 346–354.
- 417 [11] I. Carpinteiro, M. Ramil, I. Rodríguez, R. Cela, *J. Chromatogr. A* 1217 (2010) 7484–7492.
- 418 [12] P. Cabras, A. Angioni, V.L. Garau, F.M. Pirisi, V. Brandolini, *J. AOAC Int.* 81 (1998) 1185–1189.
- 419 [13] A. Angioni, M. Schirra, V.L. Garau, M. Melis, C.I.G. Tuberoso, P. Cabras, *Food Add. Contam.* 21 (2004)  
420 1065–1070.
- 421 [14] P. Viñas P, N. Campillo, N. Martínez-Castillo, M. Hernández-Córdoba, *J. Chromatogr. A* 1216 (2009)  
422 140–146.
- 423 [15] J. López-Flores, A. Molina-Díaz, M.L. Fernández-de-Córdoba, *Anal. Chim. Acta* 585 (2007) 185–191.
- 424 [16] E.D. Tsochatzis, U. Menkissoglu-Spiroudi, D.G. Karpouzias, R. Tzimou-Tsitouridou, *Anal. Bioanal. Chem.*  
425 397 (2010) 2181–2190.
- 426 [17] J. García-Reyes, A.U. Jackson, A. Molina-Díaz, R.G. Cooks, *Anal. Chem.* 81 (2009) 820–829.

- 427 [18] K. Zhang, J.W. Wong, D.G. Hayward, P. Sheladia, A.J. Krynitsky, F.J. Schenck, M.G. Webster, J.A.  
428 Amman, S.E. Ebeler, *J. Agric. Food Chem.* 57 (2009) 4019–4029.
- 429 [19] J. Schurek, L. Vaclavik, H. Hooijerink, O. Lacina, J. Poustka, M. Sharman, M. Caldow, M.W.F. Nielen, J.  
430 Hajslova, *Anal. Chem.* 80 (2008) 9597–9575.
- 431 [20] A. Sannino, L. Bolzoni, M. Baldini, *J. Chromatogr. A* 1036 (2004) 161–169.
- 432 [21] K.Y.F. Yau, H. Lee, J.C. Hall, *Biotechnol. Advances* 21 (2003) 599–637.
- 433 [22] V.S. Morozova, A.I. Levashova, S.A. Eremin, *J. Anal. Chem.* 60 (2005) 202–217.
- 434 [23] J.M. Van Emon, J.C. Chuang, K. Dill, G. Xiong, in: J.L. Tadeo (Ed.), *Analysis of Pesticides in Food and*  
435 *Environmental Samples*, CRC Press, Boca Raton, 2008, pp. 95–123.
- 436 [24] K. Landsteiner, *The Specificity of Serological Reactions*, revised ed., Dover Publications, New York,  
437 1962.
- 438 [25] G.S. Furzer, L. Veldhuis, J.C. Hall, *J. Agric. Food Chem.* 54 (2006) 688–693.
- 439 [26] J.V. Mercader, C. Suárez-Pantaleón, C. Agulló, A. Abad-Somovilla, A. Abad-Fuentes, *J. Agric. Food*  
440 *Chem.* 56 (2008) 1545–1552.
- 441 [27] J.V. Mercader, C. Suárez-Pantaleón, C. Agulló, A. Abad-Somovilla, A. Abad-Fuentes, *J. Agric. Food*  
442 *Chem.* 56 (2008) 2581–2588.
- 443 [28] J.V. Mercader, C. Suárez-Pantaleón, C. Agulló, A. Abad-Somovilla, A. Abad-Fuentes, *J. Agric. Food*  
444 *Chem.* 56 (2008) 7682–7690.
- 445 [29] F.A. Esteve-Turrillas, J. Parra, C. Agulló, A. Abad-Somovilla, A. Abad-Fuentes, J.V. Mercader, *Anal.*  
446 *Chim. Acta* 682 (2010) 93–103.
- 447 [30] J. Parra, J.V. Mercader, C. Agulló, A. Abad-Fuentes, A. Abad-Somovilla, *Tetrahedron* 67 (2011) 624–635.
- 448 [31] A. Abad, J.J. Manclús, F. Mojarrad, J.V. Mercader, M.A. Miranda, J. Primo, V. Guardiola, A. Montoya, J.  
449 *Agric. Food Chem.* 45 (1997) 3694–3702.
- 450 [32] W.L. Shelver, D.J. Smith, E.S. Berry, *J. Agric. Food Chem.* 48 (2000) 4020–4026.
- 451 [33] C.K. Jana, E. Ali, *Steroids* 64 (1999) 228–232.
- 452 [34] A. Boudina, C. Emmelin, A. Baaliouamer, O. Paissé, J.M. Chovelon, *Chemosphere* 68 (2007) 1280–1288.

453 [35] Joint FAO/WHO Meeting on Pesticide Residues (JMPR) (2008) at [http://www.fao.org/](http://www.fao.org/fileadmin/templates/agphome/documents/Pests_Pesticides/JMPR/Evaluation08/Azoxystrobin.pdf)  
454 [fileadmin/templates/agphome/documents/Pests\\_Pesticides/JMPR/Evaluation08/Azoxystrobin.pdf](http://www.fao.org/fileadmin/templates/agphome/documents/Pests_Pesticides/JMPR/Evaluation08/Azoxystrobin.pdf).  
455 Accessed 14 April 2011.

456 [36] J.M. Clough, C.R.A. Godfrey, J.J. Eshelby, GB Patent 2248615(A), 1992.

457

458

459

460

461

462

463 **Figure legends**

464 **Fig. 1.** Chemical structures of AZ and haptens used for immunization and assay development.

465 **Fig. 2.** Minimum, mean, and maximum CR values for the *Z*-isomer of AZ exhibited by the four  
466 different mAb populations.

467 **Fig. 3.** Influence of buffer pH and ionic strength on the  $A_{\max}$  and  $IC_{50}$  values of the selected assay.

468 Antibody AZo6#49 was used at  $60 \text{ ng mL}^{-1}$  and the OVA–AZb6 conjugate was used at  $1.0 \text{ }\mu\text{g mL}^{-1}$ .

469 Values are the mean of three independent experiments.

470 **Fig. 4.** Matrix effects for tomato, peach, grape, and banana juices. Food samples were diluted in buffer

471  $1/5$  (●),  $1/25$  (○),  $1/50$  (◆),  $1/250$  (Δ),  $1/500$  (■). A control was also run with no juice (□).

472

**Table 1**  
Antibody evaluation and assay selection using AZ as competitor in the i-cELISA format

mAb	Coating conjugate															
	OVA-AZa6				OVA-AZb6				OVA-AZc6				OVA-AZo6			
	[mAb] <sup>a</sup>	[C] <sup>b</sup>	slope	IC <sub>50</sub> <sup>c</sup>	[mAb]	[C]	slope	IC <sub>50</sub>	[mAb]	[C]	slope	IC <sub>50</sub>	[mAb]	[C]	slope	IC <sub>50</sub>
AZa6#11	30	1.0	-1.23	0.25	20	1.0	-1.28	0.14	1000	1.0	<sup>d</sup>		30	0.1	-1.83	0.30
AZa6#21	100	0.1	-1.59	0.87	25	1.0	-1.24	0.65	300	1.0	-1.73	1.66	30	0.1	-1.34	0.36
AZa6#210	100	0.1	-1.35	0.44	1000	1.0			1000	1.0			30	0.1	-1.11	0.42
AZa6#33	30	0.1	-1.76	0.23	30	1.0	-1.51	0.21	1000	1.0			20	0.1	-1.81	0.15
AZb6#22	1000	1.0			100	0.1	-1.34	0.35	100	1.0	-1.49	0.96	100	0.1	-1.36	0.42
AZb6#23	1000	1.0			100	0.1	-1.43	0.79	100	1.0	-1.17	1.27	100	0.1	-1.31	0.91
AZb6#24	1000	1.0			100	0.1	-1.34	0.29	100	1.0	-1.37	0.88	100	0.1	-1.40	0.36
AZb6#38	300	0.1	-1.17	6.97	100	0.1	-1.26	1.99	1000	1.0			30	0.1	-0.98	3.94
AZc6#22	100	1.0	-1.63	0.39	100	0.1	-1.81	0.47	100	0.1	-1.59	0.40	100	0.1	-1.47	0.39
AZc6#24	100	1.0	-1.42	0.53	300	0.1	-1.64	0.31	100	0.1	-1.64	0.30	100	0.1	-1.43	0.31
AZc6#27	100	0.1	-1.39	0.54	300	1.0	-1.61	1.70	100	0.1	-1.47	0.62	30	0.1	-1.13	0.46
AZc6#31	1000	1.0			1000	1.0			100	0.1	-0.95	0.68	100	0.1	-1.08	0.49
AZo6#31	100	0.1	-1.87	1.28	15	1.0	-1.57	0.20	1000	1.0			30	0.1	-1.89	0.38
AZo6#44	100	0.1	-1.95	0.45	30	1.0	-1.24	0.21	300	1.0	-1.62	0.79	50	0.1	-1.50	0.25
AZo6#45	100	0.1	-1.93	0.67	1000	1.0			30	1.0	-1.59	0.18	30	0.1	-1.40	0.45
AZo6#49	100	0.1	-1.58	0.60	30	1.0	-1.23	0.29	300	1.0	-1.57	0.95	30	0.1	-1.45	0.29

<sup>a</sup> Antibody concentration in ng mL<sup>-1</sup>.

<sup>b</sup> Coating conjugate concentration in µg mL<sup>-1</sup>.

<sup>c</sup> Values in nM.

<sup>d</sup> Empty cells indicate that no signal was obtained for this particular combination.

**Table 2**

Curve parameters of the selected assays


mAb	Coating hapten	Parameters <sup>a</sup>		
		A <sub>max</sub>	Slope	IC <sub>50</sub> (nM)
AZa6#11	AZa6	1.3 ± 0.3	-1.2 ± 0.2	0.23 ± 0.04
	AZb6	1.3 ± 0.3	-1.3 ± 0.2	0.15 ± 0.06
AZb6#24	AZb6	0.8 ± 0.2	-1.3 ± 0.2	0.33 ± 0.03
AZc6#24	AZo6	1.3 ± 0.3	-1.3 ± 0.2	0.40 ± 0.05
AZo6#49	AZo6	1.3 ± 0.3	-1.2 ± 0.3	0.29 ± 0.05
	AZb6	1.2 ± 0.2	-1.3 ± 0.1	0.23 ± 0.03

<sup>a</sup> Mean values of five independent experiments. The lower asymptotes (A<sub>min</sub>) were always below 0.06.



**Table 3**

Standard curve, conditions, and parameters of the developed assay for AZ<sup>a</sup>



Antibody	mAb AZo6#49 at $60 \text{ ng mL}^{-1}$
Conjugate	OVA-AZb6 at $1.0 \mu\text{g mL}^{-1}$
$A_{\text{max}}$	$1.0 \pm 0.1$
$A_{\text{min}}$	$0.026 \pm 0.007$
Slope	$-1.15 \pm 0.08$
$\text{IC}_{50}$ ( $\mu\text{g L}^{-1}$ )	$0.10 \pm 0.01$
LOD ( $\mu\text{g L}^{-1}$ )	$0.017 \pm 0.002$
Working range ( $\mu\text{g L}^{-1}$ )	$0.035 \pm 0.004$ to $0.46 \pm 0.05$
Assay buffer	PBST
Assay time (h)	2.5
Solvent tolerance	methanol 10%

<sup>a</sup> Values are the means of 16 independent determinations.

**Table 4**Recovery values (%) from spiked samples as determined by i-cELISA<sup>a</sup>

Sample	Spiked ( $\mu\text{g L}^{-1}$ )	Dilution factor		
		1/25	1/100	1/500
Tomato juice	2	115 ± 14	101 ± 8	--- <sup>b</sup>
	5	107 ± 5	102 ± 21	---
	10	123 ± 15	117 ± 11	98 ± 10
	50	67 ± 25	105 ± 12	97 ± 8
	200	---	104 ± 13	109 ± 16
	500	---	101 ± 11	106 ± 15
Peach juice	2	103 ± 26	127 ± 22	---
	5	105 ± 18	106 ± 20	---
	10	119 ± 16	118 ± 16	105 ± 18
	50	91 ± 12	111 ± 12	104 ± 13
	200	---	103 ± 15	106 ± 6
	500	---	88 ± 11	103 ± 7
Grape must	2	113 ± 16	113 ± 22	
	5	113 ± 13	99 ± 9	
	10	121 ± 18	112 ± 16	98 ± 16
	50	84 ± 23	105 ± 13	97 ± 17
	200		101 ± 10	99 ± 7
	500		97 ± 13	102 ± 9
Banana juice	2	131 ± 23		
	5	114 ± 18	118 ± 19	
	10	133 ± 20	119 ± 21	122 ± 3
	50	83 ± 13	107 ± 8	102 ± 7
	200		93 ± 19	104 ± 6
	500		96 ± 27	98 ± 6

<sup>a</sup> Values are the mean of six independent determinations.<sup>b</sup> Concentration out of range.

Figure 1  
[Click here to download high resolution image](#)




Figure 2  
[Click here to download high resolution image](#)




Figure 3

[Click here to download high resolution image](#)





Figure 4  
[Click here to download high resolution image](#)

