Development of a radiobiology beam line at the 18 MeV proton cyclotron facility at CNA (Seville, Spain)

<u>A. Baratto-Roldán^{1,2,*}, M. C. Battaglia², M. A. Cortés-Giraldo², M. C. Jiménez-Ramos¹, J. García-López^{1,2}, M. I. Gallardo², J. M. Espino^{1,2}</u>

(1) Centro Nacional de Aceleradores (Seville, Spain)
(2) Universidad de Sevilla (Seville, Spain)
(*) <u>abaratto@us.es</u>

2nd OMA Topical Workshop – Diagnostics for Beam and Patient Monitoring

CERN, Switzerland June 4th-5th, 2018








BIOLOGICAL EFFECT OF PROTON BEAMS

Preparation of beam lines at the National Centre of Accelerators (CNA) in Seville for radiobiology experiments with LOW ENERGY PROTON BEAMS



• Uniform RBE of 1.1 assumed in treatment planning with SOBP. Radiobiology experiments at low proton energies are necessary to understand how RBE varies with LET.

• **RBE increases** with LET.





Anna Baratto-Roldán 2nd OMA Topical Workshop

CNA PROTON FACILITIES



Pelletron 9SDH-2 Tandem Accelerator

Maximum terminal voltage 3MV


Maximum proton energy of 6 MeV


 SNICS II source: sputtering Cesium source (sputtering from solid targets)

Cyclone 18/9 Cyclotron Accelerator

- Maximum proton energy of 18 MeV
- Extracted maximum beam intensity: 80 μA ± 10%
- Used for production of PET radioisotopes








1. Beam Centering: Beam Profile Monitors (BPMs) with high beam current










- 1. Beam centering: Beam Profile Monitors (BPMs) with high beam current;
- 2. Decreasing the current (tenths of nA in the Faraday Cup FC):
 - a. Defocusing the beam with the quadrupole magnets (fast procedure);
 - b. Decreasing the pressure of the gas used for the stripping process (slow method)




Beam optics unaltered and the maximum beam intensity maintained at the same position





M. Cristina Battaglia, Anna Baratto-Roldán 2nd OMA Topical Workshop






- 1. Beam centering: Beam Profile Monitors (BPMs) with high beam current;
- 2. Decreasing the current (tenths of nA in the Faraday Cup FC);
- 3. Homogeneity: Beam scattering on heavy targets





M. Cristina Battaglia, Anna Baratto-Roldán 2nd OMA Topical Workshop





- 1. Beam centering: Beam Profile Monitors (BPMs) with high beam current;
- 2. Decreasing the current (tenths of nA in the Faraday Cup FC);
- 3. Homogeneity: Beam scattering on heavy targets;
- 4. Setup for sample irradiation.








SETUP FOR SAMPLE IRRADIATION

- ➤ Thin vacuum exit window (kapton, thickness 50 µm, diameter 44 mm, P≈10⁻⁶ mbar vacuum pressure).
- Ionization chamber (IC)
 - three parallel electrodes 7.5 μm thick;
 - > Two air gaps 6.75 mm thick;
 - ➢ V_{IC}=400 V;
 - connected to a Keithley electrometer.
- Holder with six positions designed for biological samples and also used for radiochromic films.







CHARACTERIZATION OF THE BEAM PROFILE WITH EBT3 FILMS


EBT3 Radiochromic fims:



EBT3 Films irradiated with different doses

Mainly used for photon dosimetry:

- Under ionizing radiation, the sensitive gel layer polimerizes, and the film turns blue in an amount that is proportional to the dose.
- Symmetric construction.
- Energy independence.
- High spatial resolution (25μm).
- Tissue equivalence.
- No chemical, thermal or optical development.







DOSIMETRY WITH EBT3 FILMS

EBT3 Radiochromic fims:

 $D_{IC} = F_{IC}$

 $F_{IC} = Q_{IC}$



EBT3 Films irradiated with different doses **Proton dosimetry:**


- Energy independence not necessarily true;
- Ionization chamber used to measure fluence and evaluate absorbed dose.

dE

 $\mathfrak{I}_{Lu}\overline{\mathrm{dx}}$

W

 ΔE_{IC}



M. Cristina Battaglia, Anna Baratto-Roldán 2nd OMA Topical Workshop

SIMULATION



METHODOLOGY

- A dedicated setup which provides the best conditions to irradiate biological samples with low energy protons (E = 3.8 MeV).
- A protocol for the beam optimization has been established, to obtain low current beams and homogeneous beam profiles .


RESULTS

- Outside the Bragg peak region: dose calibration for the EBT3 film under photon irradiation, can be successfully adopted for low energy protons, as verified through the IC fluence measurements.
- Degrading the beam energy: effects of saturation in the EBT3 film response (LET increases). This effect is more remarkable for doses higher than 10 Gy.
- The protocol for the beam optimization and the dosimetry studies are necessary for the irradiation of cell cultures. First measurements were performed with breast cancer cells.






EXTENSION TO HIGHER ENERGIES: CYCLOTRON FACILITY




CYCLOTRON FACILITY: CHALLENGES

- Experimental beam line common to different research groups and devoted to interdisciplinary research activities.
- Beam time limited by normal workflow of PET radioisotopes production.
- Limited space in the experimental room.
- External beam in-air.







- 1. Beam centering:
 - a. High beam current and focalized beam for centring procedure;
 - b. Beam position checked with scintillating foils attached to the Faraday Cup and placed at the sample position.









- 1. Beam centering;
- 2. Beam current:
 - a. Ion source current set to the minimum \sim 5 mA;
 - b. Completely defocused beam (detuned quadrupoles);
 - c. Current in the Faraday Cup of the order of hundreds of pA;
 - d. Approximately 80% of the beam intensity is lost after second collimation.




- 1. Beam centering;
- 2. Beam current;
- 3. Beam shape:
 - a. Two collimators with diameter 1.5 cm at a distance of 1 m;
 - b. Approximately elliptical (a \approx 15 mm, b \approx 12 mm).








b. Variable distance between exit window and sample.



- 1. Beam centering;
- 2. Beam current;
- 3. Beam shape;
- 4. Homogeneity;
- 5. Setup for sample irradiation.










DOSIMETRY WITH EBT3 FILMS: PHOTON CALIBRATION


PHOTON CALIBRATION @HOSPITAL: "Hospital Universitario Virgen Macarena"




- Films irradiated with a 6 MV linac under reference conditions (0 20 Gy);
- Calibration curves of the form: $D = a + b \cdot \text{netOD} + c \cdot \text{netOD}^d$



FULL-BEAM PROFILES



FULL-BEAM PROFILES



2nd OMA Topical Workshop


CHARACTERIZATION OF THE BEAM PROFILE


1. Decreasing the current: minimum beam intensity settable (pA) + defocusing;

2. Homogeneity: Beam scattering on heavy targets placed after the exit window in air.

Anna Baratto-Roldán

2nd OMA Topical Workshop






of the sample (σ): X : 1.99 ± 0.04 Gy Y : 2.04 ± 0.04 Gy Max deviation from mean value: Χ:2.8 σ Υ:2.7 σ


CHARACTERIZATION OF THE BEAM PROFILE

1. Decreasing the current: minimum beam intensity settable (pA) + defocusing;

2. Homogeneity: Beam scattering on heavy targets placed after the exit window in air.

2nd OMA Topical Workshop






Max deviation from mean value:

Χ:2.3 σ

Υ:2.2 σ

DOSIMETRY WITH EBT3 FILMS: PROTON CALIBRATION

Proton dose evaluation with the ionization chamber coupled with EBT3 films






Anna Baratto-Roldán 2nd OMA Topical Workshop


DOSIMETRY WITH EBT3 FILMS: PROTON CALIBRATION

Proton dose evaluation with the ionization chamber coupled with EBT3 films






MONTE CARLO SIMULATIONS: SRIM



2nd OMA Topical Workshop






MONTE CARLO SIMULATIONS: COMPARISON SRIM – GEANT4

125 μm PET + 150 μm W + 49.9 cm Air + IC + EBT3

IONIZATION CHAMBER






Anna Baratto-Roldán 2nd OMA Topical Workshop

MONTE CARLO SIMULATIONS: COMPARISON SRIM – GEANT4

125 μm PET + 150 μm W + 49.9 cm Air + IC + EBT3


ENERGY DEPOSITED WITHIN ACTIVE VOLUME OF THE EBT3 FILM (LUCITE)



MONTE CARLO SIMULATIONS: COMPARISON SRIM – GEANT4

125 μ m PET + 150 μ m W + 49.9 cm Air + IC + EBT3

FLUENCE SPECTRA AND AVERAGE LET IN EBT3 FILM (LUCITE)





Preparation of beam line: Next steps

- Accurate measurement of beam energy at the position of IC and EBT3 films, to be compared with SRIM and Geant4 simulations;
- Understand source of the differences between Geant4 and SRIM simulations;
- Complete Geant4 simulation of the beam characteristics and compare simulation results with full beam profile mesurements;
- Improve proton calibration with EBT3 films.
- Improve beam profile homogeneity (insertion of scattering foils in vacuum, change collimation system...)



- Thanks to Dr. M. Cristina Battaglia, who generously gave me material for this talk.
- This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Marie Sklodowska-Curie grant agreement No 675265, OMA – Optimization of Medical Accelerators, and from the Spanish Ministry of Economy and Competitiveness under grant No FPA2016-77689-C2-1-R. The Monte Carlo simulations were carried out at the FIS-ATOM cluster hosted at CICA (Seville, Spain).





THANK YOU FOR YOUR ATTENTION!



SPARE SLIDES



Goal of the project and CNA facilities;

Previous experiments at Tandem;

Extension to higher energies: the Cyclotron;


Results;


Conclusions.



CHARACTERIZATION OF THE BEAM PROFILE

Preliminary analysis of radiochromic films to identify the best Au scatterer:





Raw response for three values of irradiation with a gold target of 2.0 mg/cm²: curved profile, not suitable for the measurements





M. Cristina Battaglia, Anna Baratto-Roldán 2nd OMA Topical Workshop



CHARACTERIZATION OF THE BEAM PROFILE

Preliminary analysis of radiochromic films to identify the best Au scatterer:






Raw response for three values of irradiation with a gold target of 5.4 mg/cm²: flat profile, **suitable** for the measurements



M. Cristina Battaglia, Anna Baratto-Roldán 2nd OMA Topical Workshop





- 1. Beam centering;
- 2. Beam current;
- 3. Beam shape;
- 4. Sample centering:
 - a. Samples placed on a moving table with millimetric precision;
 - b. Samples centred at the beam position.



SETUP FOR FULL BEAM PROFILES MEASUREMENT

