

Metadata of the chapter that will be visualized online

Chapter Title	Studies of the Secretory Machinery Dynamics by Total Internal Reflection Fluorescence Microscopy in Bovine Adrenal Chromaffin Cells
Copyright Year	2018
Copyright Holder	Springer Science+Business Media, LLC, part of Springer Nature
Author	Family Name Villanueva Particle Given Name José Suffix Division Instituto de Neurociencias Organization Centro Mixto CSIC-Universidad Miguel Hernández Address Alicante, Spain
Author	Family Name Gimenez-Molina Particle Given Name Yolanda Suffix Division Instituto de Neurociencias Organization Centro Mixto CSIC-Universidad Miguel Hernández Address Alicante, Spain
Corresponding Author	Family Name Gutiérrez Particle Given Name Luis M. Suffix Division Instituto de Neurociencias Organization Centro Mixto CSIC-Universidad Miguel Hernández Address Alicante, Spain Email luisguti@umh.es
Abstract	Cultured bovine chromaffin cells have been tested as a successful neuroendocrine model to study the secretory process. Changes in the dynamics of the secretory vesicles and the exocytotic machinery microdomains could be studied in control and stimulated conditions using appropriate molecular tools such as fluorescent SNARE protein expression or fluorochrome vesicular labeling in these neuroendocrine cells. Since most of these changes occur in or near the plasma membrane, the use of the total internal reflection fluorescent microscopy (TIRFM) and the implement of particle motion analysis could be essential tools to study the structural and dynamic changes of secretory machinery related with its function in this exocytotic cell model.
Keywords (separated by ‘-’)	TIRFM - Evanescent field - Exocytotic events - SNARE proteins - MSD - Diffusion coefficient - Chromaffin granules - Chromaffin cells

Studies of the Secretory Machinery Dynamics by Total Internal Reflection Fluorescence Microscopy in Bovine Adrenal Chromaffin Cells 2 3 4

José Villanueva, Yolanda Gimenez-Molina, and Luis M. Gutiérrez 5 AU1

Abstract 6

Cultured bovine chromaffin cells have been tested as a successful neuroendocrine model to study the secretory process. Changes in the dynamics of the secretory vesicles and the exocytotic machinery microdomains could be studied in control and stimulated conditions using appropriate molecular tools such as fluorescent SNARE protein expression or fluorochrome vesicular labeling in these neuroendocrine cells. Since most of these changes occur in or near the plasma membrane, the use of the total internal reflection fluorescent microscopy (TIRFM) and the implement of particle motion analysis could be essential tools to study the structural and dynamic changes of secretory machinery related with its function in this exocytotic cell model.

Key words TIRFM, Evanescent field, Exocytotic events, SNARE proteins, MSD, Diffusion coefficient, Chromaffin granules, Chromaffin cells 15
16

1 Introduction 17

Neurosecretion involves the transport of specialized vesicles to the immediate vicinity of active sites and they regulated fusion. In this sense, the initial transport of the vesicles is governed by the activity of the cell cytoskeleton [1] whereas SNARE proteins appear to be essential for the final exocytotic fusion steps [2]. Since these final stages of neurosecretion occur nearby the plasma membrane, the study of these structural and dynamic changes requires microscopic techniques with a good signal/noise ratio in the z -axis perpendicular to the plasma membrane plane. Confocal techniques have a relatively good z -axis resolution but the laser excitation field extends up to ~ 600 nm generating an elevated z -axis background [3]. On the other hand, in TIRFM techniques the excitation field, named evanescent wave, is thinner extending around ~ 200 nm in the interior of the cells attached to the coverslip surface (Fig. 1) (see Note 1) [4], and in consequence only fluorophore molecules 32




Fig. 1 Total internal reflection microscopy. When laser light (IN) strikes the interface between two optical media of different refractive indices (i.e., objective/specimen), at an angle greater than the critical angle undergoes total reflection. Beyond the angle of total reflection, the electromagnetic field (termed *evanescent wave*) still extends a few hundred nanometers into the z direction but decreasing exponentially with the distance, therefore exciting fluorescent molecules within this evanescent field (see **Note 4**)

located in the surface of the specimen are excited to emit. In other words, TIRFM makes possible to obtain a high contrast observation of the secretory machinery behavior in or near the plasma membrane having a good signal/noise ratio [5, 6].

The employment of TIRFM techniques for analyzing dynamic changes in SNARE fluorescent constructs (i.e., GFP-SNAP-25) expressed in chromaffin cells [7] has been demonstrated previously to be a useful tool for studying both the fusion complex assembly and the vesicle fusion kinetics [4]. In the same way, the use of fluorescent markers for acidic chromaffin granules has allowed to design and perform dynamic TIRFM assays to analyze both the mobility and fusion parameters changes in response to different secretory conditions [4, 8, 9].

This work describes the methodological approach in order to carry out these studies in bovine chromaffin cells in monolayer cultures.

2 Materials

2.1 Basal and Depolarizing Solutions

All solutions have to be prepared using ultrapure water (purifying deionized water, sensitivity of 18 MΩ-cm at 25 °C and analytical grade reagents). Prepare and store all reagents at 4 °C (unless indicated otherwise).

1. Krebs/HEPES (K/H) basal solution (134 mM NaCl, 4.7 mM KCl, 1.2 mM KH₂PO₄, 1.2 mM MgCl₂, 2.5 mM CaCl₂,

	11 mM glucose, 0.56 mM ascorbic acid, and 15 mM HEPES, pH 7.4).	56
		57
	2. Krebs/HEPES (High K) depolarizing solution (80 mM NaCl, 59 mM KCl, 1.2 mM KH_2PO_4 , 1.2 mM MgCl_2 , 2.5 mM CaCl_2 , 11 mM glucose, 0.56 mM ascorbic acid, and 15 mM HEPES, pH 7.4).	58
		59
		60
		61
		62
2.2 Regulation for the Use of Bovine Adrenal Glands	Bovine adrenal glands were obtained from an industrial slaughterhouse that is subject to strict regulations issued by the Spanish Ministries of Agriculture, Industry and Health, and in accordance with European Community guidelines.	63
		64
		65
		66
	All the protocols described here were approved by the “Organo Evaluador de Proyecto” at the University Miguel Hernández, the office in charge of overseeing the ethical issues associated with animal care and experimentation at our investigation institute.	67
		68
		69
		70
		71
2.3 Chromaffin Cell Reagents, Solutions, and Fungibles	1. Locke buffer $1\times$ (154 mM NaCl, 5.58 mM KCl, 3.6 mM NaHCO_3 , 5.6 mM glucose, 5 mM HEPES, sterilized, adjusted to pH 7.4, and stored at 4 °C).	72
		73
		74
	2. Type A Collagenase (store at 4 °C).	75
		76
	3. Bovine serum albumin (BSA) (store at 4 °C).	76
		77
	4. Type A Collagenase Solution (Cold Locke buffer $1\times$, 0.25% Type A Collagenase, 0.5 % BSA) should be prepared before use.	77
		78
	5. Percoll continuous gradient (store at 4 °C).	79
		80
	6. Diluted Percoll Continuous Gradient (mix fresh 36 mL of Percoll Continuous Gradient with 4 mL of Locke buffer $1\times$, maintaining at 37 °C until use).	81
		82
	7. Dulbecco’s modified Eagle’s medium (DMEM) (store at 4 °C).	83
		84
	8. Supplemented DMEM (DMEM supplemented with 10% fetal calf serum, 10 μM cytosine arabinoside, 10 μM 5-fluoro-2’-deoxyuridine, 50 IU/mL penicillin, and 50 $\mu\text{g}/\text{mL}$ streptomycin, stored at 4 °C).	85
		86
		87
		88
	9. 35 mm Petri dishes.	89
		90
2.4 Fluorescent Protein Expression	The use of expression vectors is a very useful tool for dynamic studies using fluorescence microscopy. In this way, we used the GFP-SNAP-25 construct like an example to analyze SNARE protein behavior. This construct was obtained by cloning the cDNA corresponding to the SNAP-25a isoform [10] into the XhoI and BamHI sites of pEGFP-C3 expression vector (Clontech, Palo Alto, CA). This construct expresses the protein SNAP-25 fused in-frame at the C-terminus to EGFP.	91
		92
		93
		94
		95
		96
		97
		98
		99

2.5 Fluorescent Beads, Fluorochrome Dyes, and Working Solutions	1. 100 nm Fluorescent beads.	100
	2. Lyso Red Stock Solution: LysoTracker Red DND-99 (Invitrogen). Resuspend at 1 mM in DMSO.	101
	3. Acridine Orange Stock Solution: Acridine Orange.	102
	4. Lyso Green Stock Solution: LysoTracker Green DND-26 (Invitrogen): Resuspend at 1 mM in DMSO.	103
	5. Prepare Lyso Red or Lyso Green Working Solutions (1 μ M) by adding 1 μ L from each stock solution to a 1 mL of Krebs/HEPES (K/H) basal solution.	104
	6. Prepare Acridine Orange Working Solution (2 μ M) by adding 2 μ L from Acridine Orange Stock Solution to a 1 mL of Krebs/HEPES (K/H) basal solution.	105
2.6 DNA Electroporation Kit	Amaxa basic nucleofector kit for primary mammalian neuronal cells (Amaxa GmbH).	106
2.7 TIRFM Microscopy	1. A through-the-lens TIRFM system using an Olympus IX-71 inverted microscope with a 100 \times PlanApo oil immersion TIRFM objective with 1.45 N.A. (numeric aperture).	107
	2. Epifluorescence and laser illumination (488 nm argon ion 40 mW or 543 nm He/Ne 10 mW) using an Olympus TIRFM IX2-REAEVA combiner system that allows changing of the angle of laser incidence.	108
	3. Fluorescence emission is split using an Optosplit II system (Cairn Research Ltd.) equipped with GFP and rhodamine filter sets.	109
	4. Acquire separated images simultaneously using an Electron Multiplier CCD cooled camera (Hamamatsu), controlled by the Imaging Software Wasabi v.1.5 (Hamamatsu Photonics) in an IBM-compatible PC.	110
2.8 MSD Algorithm and Analysis Software	1. Use the public domain ImageJ program with plug-ins to analyze the fluorescent images.	111
	2. After thresholding of the images, determine particle centroid in time-lapse studies using a multi-tracker plug-in as described earlier [3, 11, 12].	112
	3. Transform the x - y coordinates for centroids of these particles to Igor Pro program (WaveMetrics Inc.), and use specialized macros to calculate the total lateral displacement and the mean square displacement (MSD) for any given time interval by using Eqs. (1) and (2) (see Notes 1 and 2).	113

3 Methods

142

3.1 Isolation, Culture, and Transfection of Bovine Chromaffin Cells

1. Prepare bovine chromaffin cells by collagenase digestion and 143
separate from debris and blood cells using a centrifugation in 144
Percoll gradients as described before [13]. 145
2. Store adrenal glands in cold (4 °C) Locke buffer. 146
3. Trim each gland of fat and perfuse, using a syringe, via adrenal 147
vein with 3–5 mL of 37 °C warmed Type A Collagenase Solu- 148
tion over a period of 15 min while maintained at 37 °C. Repeat 149
up to three times. 150
4. Dissect each gland using scissors and mechanically scrape and 151
disaggregate the exposed medulla to obtain a cell homogenate. 152
5. Resuspend the resulting solution with warmed (37 °C) Locke 153
buffer and filter successively using first a sieve of 217 microme- 154
ter diameter and a second one of 82 µm diameter. 155
6. Add the resulting cell suspension to a Locke buffer to a final 156
volume of 200 mL and centrifuge ($100 \times g$ for 5 min) in Locke 157
buffer at room temperature to wash out collagenase. 158
7. Resuspend cell pellets using Locke buffer up to a final volume 159
of 42 mL and mix with 38 mL of diluted Percoll continuous 160
gradient. 161
8. Centrifuge for 22 min at 37 °C in a Beckman JA 25.50 rotor at 162
 $21,000 \times g$. 163
9. Collect cells equilibrated between the densities of 1.045 and 164
1.075 g/mL. This gradient fraction contains chromaffin cells 165
without blood cells and debris. 166
10. Wash the isolated chromaffin cells twice by centrifugation 167
($100 \times g$ for 5 min) in Locke buffer at room temperature. 168
11. Harvest isolated chromaffin cells and resuspend with 10 mL of 169
Dulbecco's modified Eagle's medium (DMEM). 170
12. Count cells in a Neubauer chamber to estimate cell number 171
and concentration. 172
13. To obtain culture cells, plate aliquots containing 750,000 173
isolated chromaffin cells in 35 mm Petri dishes as monolayer 174
cultures with supplemented DMEM ($150,000 \text{ cells/cm}^2$), and 175
transferred to an incubator 37 °C 5% CO₂. These cells could be 176
used for experimentation during a week. 177
14. Transfect aliquots containing five million isolated chromaffin 178
cells with 1–3 µg of the appropriate DNA expression plasmid 179
(GFP-SNAP-25) using the Amaxa basic nucleofector kit for 180
primary mammalian neuronal cells according to the manufac- 181
turer's instructions. 182

	15. Plate transfected cells in 35 mm Petri dishes as monolayer cultures with supplemented DMEM (150,000 cells/cm ²), and maintained in a temperature incubator as described before (<i>see step 4</i>). Use the transfected cells between the second and fourth days after plating.	183 184 185 186 187 188
3.2 Fluorochrome Vesicle Labeling	1. Replace the culture media in each 35 mm Petri dish with Krebs/HEPES (K/H) basal solution.	189 190
	2. To stain a population of acidic vesicles for experiments involving the motion of granules incubate cells for 20 min at room temperature and darkness with either Lyso Red or Lyso Green Working solution (1 μM) [14].	191 192 193 194
	3. For experiments aiming to measure both vesicular motion and fusion incubate cells for 20 min at room temperature and darkness with Acridine Orange Working Solution (2 μM) [4].	195 196 197
	4. Wash cells extensively with Krebs/HEPES (K/H) basal solution and mount to use them for TIRFM experiments within the next 2 h.	198 199 200 201
3.3 TIRFM Laser Calibration and Alignment	1. Calibrate TIRFM using 100 nm fluorescent beads adhered to the coverslip.	202 203
	2. Determine the fluorescence intensities of the beads at different vertical planes with step sizes of 100 nm using the motorized system mounted on the microscope. Obtain images for both epifluorescence and TIRFM.	204 205 206 207
	3. Estimate the depth of penetration for the evanescent field to ~200 nm ($d = 180 \pm 16$ nm) (<i>see Note 5</i>) to permit visualization of the static beads adhered to the coverslip. Beads in suspension undergoing random movement should be infrequently seen in TIRFM, and the vast majority visualized by epifluorescence.	208 209 210 211 212 213 214
3.4 TIRFM Dynamic Assays	1. For vesicle motion assays, label cells with red acridine orange fluorescence or by LysoTracker red fluorescence (<i>see Subheading 3.2</i>). Visualize the cells using a laser excitation 543 nm He/Ne 10 mW [9].	215 216 217 218
	2. For SNARE position, use cells expressing GFP-SNAP-25 (<i>see Subheading 3.1</i>). Visualize cells using laser excitation 488 nm argon ion 40 mW [4, 9].	219 220 221
	3. Wash cells with Krebs/HEPES (K/H) basal solution for control conditions or treat with different substances for a variety of experimental conditions. After extensive washing, observe the cells under TIRFM using the appropriate laser excitation (543 nm He/Ne 10 mW for labeled vesicles or 488 nm argon ion 40 mW for GFP-SNAP25). Acquire fluorescence images at 1-s intervals during 20 s.	222 223 224 225 226 227 228




Fig. 2 TIRFM dynamic assays. Time-lapse original fluorescent images were processed and imported in ImageJ program to select particles and measure centroid XY displacement (*see* Subheading 2.4). XY coordinates for centroids were imported and processed with specialized macros in Igor Pro program to calculate the mean square displacement (MSD) and the diffusion parameters for any given time interval [15] (*see* **Note 6**)

4. For the motion analysis, threshold fluorescent images to deter- 229
mined particle centroid in time-lapse studies using a plug-in in 230
Igor Pro program as described earlier. Analyze trajectories to 231
calculate average speeds, in each assay conditions. 232
5. Use specialized macros in Igor Pro program *xy* coordinates for 233
centroids to calculate the total lateral displacement and the 234
mean square displacement (MSD). 235
6. Export MSD data to Prism Graph Pad v.4.01 software to 236
represent the MSD curves (Fig. 2) [15] (*see* **Note 6**). In each 237
condition of assay, express the averaged MSD data from a 238
number (*n*) of individual cells from at least three different 239
cultures as a mean \pm S.E.M. 240
7. Calculate the diffusion coefficient *D* (*see* **Note 2**, Eq. (2)) in 241
some assays assuming that movement was governed by a single 242
coefficient of diffusion and its value could be derived from the 243
fitted slope of the lineal regression for averaged MSD data: 244
 $D = \text{slope}/4$ [4, 16]. 245
8. In other assays, fit MSD-averaged data to a nonlinear regres- 246
sion using a hyperbolic equation to determine whether dis- 247
placement of the particles could be restricted to an area 248
(Fig. 2) [8]. 249
9. Establish the statistical significance of changes in experimental 250
data between different assay conditions using the Student's *t*- 251
test for paired samples or the two-way ANOVA test (behaviors 252
were considered significantly different when $P < 0.05$). 253

3.5 TIRFM Fusion Assays

1. For vesicle fusion assays, using acridine orange to label granules (see Subheading 3.2). Assess the granules by the red acridine orange fluorescence in mature acidic vesicles whereas their fusion is followed by the green flashes produced after granule matrix neutralization during exocytosis according to procedures described previously [4] or by the total or partial disappearance of red acridine orange fluorescence after exocytosis (Fig. 3).
2. Stimulate cells for 1 min using cell perfusion with Krebs/HEPES (High K) depolarizing solution.
3. In these experiments, keep the laser intensity low (2–4% of the maximal intensity with the 488 nm argon ion 40 mW laser) to prevent light-induced fusion [17]. No fusion should be detected in the absence of cell stimulation.
4. For the analysis of acridine orange-labeled vesicle fusion take images in the green channel at 20-ms intervals. Subject fusion flashes to maximal intensity determination and the data transferred to Igor Pro.
5. Analyze fusion events using software developed for amperometric detection of exocytotic events (Quanta analysis [18]). Obtain the kinetic parameters such as the time at the half-height amplitude ($t_{1/2}$) for hundreds of fusion events and represent distributions using Prism Graph Pad v.4.01 software.




Fig. 3 TIRFM fusion assays. Using a perfusion system, cells were stimulated for 1 min using a depolarizing solution with 59 mM high potassium. Green channel images showed transitory fusion flashes and red channel images showed a disappearing of previous emission. Both registers were subjected to intensity determination and transferred to Igor Pro. Fusion events were analyzed using software developed for amperometric detection of exocytotic events (Quanta analysis [18]). Kinetic parameters such as the time at the half-height amplitude ($t_{1/2}$) were obtained for hundreds of fusion events in each condition assayed

Average the fusion peak shape aged for individual cells and
examine for statistical variations (*see* Subheading 3.5, **step 9**).
Perform all the functional experiments at 21–22 °C.

4 Notes

- Equation (1)

$$\text{MSD}(n\delta t) = \frac{1}{(N-n)} \sum_{j=1}^{N-n} \left\{ [x(j\delta t + n\delta t) - x(j\delta t)]^2 + [y(j\delta t + n\delta t) - y(j\delta t)]^2 \right\} \quad (1)$$

where δt are the intervals in which images were taken. The
coordinates of a centroid were $x(j\delta t)$ and $y(j\delta t)$ for one
image and $x(j\delta t + n\delta t)$ and $y(j\delta t + n\delta t)$ for the other image.
The displacement during the interval $n\delta t$ can be measured for
($N - n$) intervals in each temporal sequence. N is the total
number of images acquired.

- Equation (2)

$$\text{MSD}(\Delta t) = r_c^2 \left[1 - A_1 \exp\left(\frac{-4A_2 D \Delta t}{r_c^2}\right) \right] \quad (2)$$

The first two terms in an infinite series are represented by this
equation. A_1 and A_2 are constants (0.99 and 0.85, respec-
tively), D is the coefficient of confined diffusion, and r_c is the
radius of the theoretical circular cage. MSD could be used to
calculate r_c when $\Delta t \rightarrow \infty$ and the data reach an asymptotic
value.

- Total reflection in the interface between two optical media with
different refractive indexes could be described by the following
Eq. (3):

$$(\theta_c) = \sin^{-1}(n(2)/n(1)) \quad (3)$$

where $n(1) > n(2)$, $n(1)$ is the refractive index of the objective,
 $n(2)$ is the refractive index of the specimen, and $\theta(c)$ is the
critical angle.

- The evanescent wave created by total reflection extends in the
 z -axis but it decreases exponentially with the distance following
this Eq. (4):

$$I = I_{(0)} e^{-z/d} \quad \text{or} \quad z = -d \ln(I/I_{(0)}) \quad (4)$$

- I is the maximum fluorescence of a particle, $I_{(0)}$ is 255 (maxi-
mum of fluorescence for 8-bit images), z is the distance from
the interphase of total reflection, and d is the constant of

exponential decay of the evanescent wave that depends on the specific properties of the TIRFM equipment (i.e., numeric aperture of the objective, etc.), by Eq. (5):

$$d = \frac{\lambda}{4\pi} (n_2^2 \sin^2 \theta - n_1^2)^{1/2} \quad (5)$$

In our system, we calculated $d = 180 \pm 16$ nm and consequently we can excite and see particles nearby plasmatic membrane until approximately 200 nm inside [4].

6. TIRFM dynamic assays can be performed using either fluorescence SNARE protein constructs or fluorescent vesicular markers. The protocol for these studies was very similar (*see* Subheading 3.4): time series TIRFM images were acquired at 1-s intervals during 20 s. The images were stored in a 16-bit HIS format, and exported in 8-bit TIFF multi-file format. Multi-files were then imported to ImageJ program and converted to an 8-bit TIFF stack time-lapse file. This file was processed using threshold option to assess the position of individual particles. The degree of threshold was specific for each type of fluorescent dye. The measurement settings were selected to find the centroid XY position for every particle studied. Black-and-white resulting images were processed with a multitracker plug-in (<https://imagej.nih.gov/ij/plugins/multitracker.html>). The minimum and maximum size parameters for particle selection and tracking were selected in base to the reported sizes found in bibliography [4, 16, 19]. Individual XY particle displacement data were copied to exported to Igor Pro datasheets. Then data were processed with the MSD macro based on MSD algorithm. Finally, MSD data were copied to export to Prism Graph Pad software v.4.01 to represent the MSD curves (Fig. 2).

Acknowledgments

This study was supported by grants from the Spanish Ministerio de Economía y Competitividad (BFU2015-63684-P, MINECO, FEDER, UE) to LMG.

References

1. Trifaro JM, Bader MF, Doucet JP (1985) Chromaffin cell cytoskeleton: its possible role in secretion. *Can J Biochem Cell Biol* 63:661–679
2. Weber T, Zemelman BV, Mcnew JA et al (1998) SNAREpins: minimal machinery for membrane fusion. *Cell* 92:759–772
3. Giner D, Neco P, Frances MM et al (2005) Real-time dynamics of the F-actin cytoskeleton during secretion from chromaffin cells. *J Cell Sci* 118:2871–2880
4. Lopez I, Ortiz JA, Villanueva J et al (2009) Vesicle motion and fusion are altered in chromaffin cells with increased SNARE cluster dynamics. *Traffic* 10:172–185

- 362 5. Oheim M, Loerke D, Stuhmer W et al (1998) 396
363 The last few milliseconds in the life of a secre- 397
364 tory granule. Docking, dynamics and fusion 398
365 visualized by total internal reflection fluores- 399
366 cence microscopy (TIRFM). *Eur Biophys J* 400
367 27:83–98 401
- 368 6. Steyer JA, Horstmann H, Almers W (1997) 402
369 Transport, docking and exocytosis of single 403
370 secretory granules in live chromaffin cells. 404
371 *Nature* 388:474–478 405
406
- 372 7. Criado M, Gil A, Viniegra S et al (1999) A 407
373 single amino acid near the C terminus of the 408
374 synaptosome-associated protein of 25 kDa 409
375 (SNAP-25) is essential for exocytosis in chro- 410
376 maffin cells. *Proc Natl Acad Sci U S A* 411
377 96:7256–7261 412
- 378 8. Villanueva J, Torres V, Torregrosa-Hetland CJ 413
379 et al (2012) F-actin-myosin II inhibitors affect 414
380 Chromaffin granule plasma membrane distance 415
381 and fusion kinetics by retraction of the cyto- 416
382 skeletal cortex. *J Mol Neurosci* 48:328–338 417
418
- 383 9. Garcia-Martinez V, Villanueva J, Torregrosa- 419
384 Hetland CJ et al (2013) Lipid metabolites 420
385 enhance secretion acting on SNARE microdo- 421
386 mains and altering the extent and kinetics of 422
387 single release events in bovine adrenal chromaf- 423
388 fin cells. *PLoS One* 8:e75845 424
425
- 389 10. Bark IC, Wilson MC (1994) Human cDNA 426
390 clones encoding two different isoforms of the 427
391 nerve terminal protein SNAP-25. *Gene* 428
392 139:291–292 429
- 393 11. Neco P, Giner D, Del Mar FM et al (2003) 426
394 Differential participation of actin- and tubulin- 427
395 based vesicle transport systems during 428
429 secretion in bovine chromaffin cells. *Eur J Neu-
ros* 18:733–742
12. Neco P, Giner D, Viniegra S et al (2004) New 398
roles of myosin II during vesicle transport and 399
fusion in chromaffin cells. *J Biol Chem* 400
279:27450–27457 401
13. Almazan G, Aunis D, Garcia AG et al (1984) 402
Effects of collagenase on the release of 403
[3H]-noradrenaline from bovine cultured 404
adrenal chromaffin cells. *Br J Pharmacol* 405
81:599–610 406
14. Becherer U, Moser T, Stuhmer W et al (2003) 407
Calcium regulates exocytosis at the level of 408
single vesicles. *Nat Neurosci* 6:846–853 409
15. Qian H, Sheetz MP, Elson EL (1991) Single 410
particle tracking. Analysis of diffusion and flow 411
in two-dimensional systems. *Biophys J* 412
60:910–921 413
16. Giner D, Lopez I, Villanueva J et al (2007) 414
Vesicle movements are governed by the size 415
and dynamics of F-actin cytoskeletal structures 416
in bovine chromaffin cells. *Neuroscience* 417
146:659–669 418
17. Henkel AW, Upmann I, Bartl CR et al (2006) 419
Light-induced exocytosis in cell development 420
and differentiation. *J Cell Biochem* 421
97:1393–1406 422
18. Mosharov EV, Sulzer D (2005) Analysis of 423
exocytotic events recorded by amperometry. 424
Nat Methods 2:651–658 425
19. Steyer JA, Almers W (1999) Tracking single 426
secretory granules in live chromaffin cells by 427
evanescent-field fluorescence microscopy. *Bio-
phys J* 76:2262–2271 428
429

Author Queries

Chapter No.: 25	442602_1_En
-----------------	-------------

Query Refs.	Details Required	Author's response
AU1	Please check and confirm if the affiliation is presented correctly.	
AU2	Please check and provide appropriate section cross-reference.	
AU3	Please check if the presentation of list is okay as typeset in Notes section.	
AU4	Please provide in-text citation for “Note 3” and also cite notes in sequential order.	
AU5	Please check if the presentation of equation is okay throughout the chapter.	
AU6	Please check the equation 5.	

Uncorrected Proof